

USING THE MICROSCOPE WITH CHROMEBOOK

Using the Southern Science Supply microscopes with Chromebooks is easy!

Most of our microscopes require a USB port for power. Many schools prefer to simply plug the microscope into the USB port and use it as an external camera on the Chromebook. Saved images can be exported to Google Draw, Google Share, etc.

IPEVO Visualizer is a robust app that works well and should be used when using a Chromebook. You can get the free app in the Chrome store, and with the handy tools such as Zoom, Rotate, Mirror, Snapshot and more, you're sure to enjoy your USB camera experience!

To learn more about the app, or to download iPEVO Visualizer follow the link:

https://www.ipevo.com/ipevo_visualizer

***NOTE:** The information provided above has not been verified by Southern Science Supply and is provided as information only. Southern Science Supply does not provide technical support for Apps.

