



SOUTHERN
SCIENCE
SUPPLY

www.southernsciencsupply.com

The MicroSight Microscope

Which scope is best for you?

First, choose the desired level of magnification that best suits your needs. **Then** select the resolution you prefer.



Magnification Options

10 – 200x...Variable focus:
best for everyday exploration

500x...Fixed focus:
for advanced learners, viewing
prepared slides and water samples

Resolution Options

1.3MP: standard resolution

5MP: high resolution,
preferred when projecting
images, cropping, printing
and making videos

5 Megapixel / 10-200x / Exclusive Software Specs

- ▶ Supported by all MAC and Windows computers, Compatible with Chromebook
- ▶ Five megapixel of pure viewing power at maximum 2592 x 1944 resolutions
- ▶ High speed USB 2.0 interface
- ▶ Wide ranges of adjustable magnification from 10X to 200X
- ▶ Viewing and capture the magnified image of your specimens right on your computer screen with VGA (640x480), XGA (1024x768) and XGA (2048x1536) resolution
- ▶ Exclusive Software with capturing picture, recording video, measuring, freezing image, digital zooming in/out, comparing
- ▶ Active pixel technology for sharp image and accurate color reproduction
- ▶ The build-in 6 White-light LED and adjustable illumination ensure the magnified images are clear and bright
- ▶ Auto exposure, auto gain control and auto white balance for optimum image
- ▶ Software controllable parameters
 - Resolution • Color • Brightness • Saturation
- ▶ High quality lens for true to life undistorted color image
- ▶ Bus powered with low power consumption - ideal for portable equipment



SOUTHERN
SCIENCE
SUPPLY

www.southernsciencesupply.com

The MicroSight Microscope